DESCRIZIONE: BENCH 4 WAY STANDOUT

APPLICATION:

Bench mounting

DETAILS:

- Test pressure: 9 Bars/air
- Max working pressure:

7 Bars

SURFACE PROTECTION:

 Epoxy powder coating RAL 7035

MATERIALS:

- Brass
- Handle in PP

HEADWORK SEALING:

CERAMIC AND Nitrile


CONNECTION TAP:

Male thread G1/2"


STANDARD AND NORMS:

- EN13792:2000
- DIN12898
- DIN12918
- ISO228/1

OPTIONS:With chemical nickel plating 17 microns inside



Codice disegno: 12013.2MDS



WEIGHT:

• 1059 gr.

Nuova FAR Rubinetterie

Disegnato da	Verificato da	
Davide Riva	Michele Albertone	
Approvato da	Data di emissione: 09/11/2011	
Michele Albertone	Peso:	Scala:1:2,5

Il presente disegno a termine di legge non può essere riprodotto o reso pubblico senza la nostra autorizzazione scritta.

DISEGNATO CON CAD NON MODIFICABILE A MANO

