

DESCRIZIONE: BENCH MOUNTED PUSH LEVER EXTRACTABLE HANDSHOWER

Codice disegno: 11090.0MDS

**APPLICATION:**

- Bench mounting

**DETAILS:**

- Waterflow around 13l/min at 3 Bars pressure
- Test pressure: 9 Bars/air
- Max working pressure:

10 Bars

**SURFACE PROTECTION:**

- Epoxy powder coating

RAL 7035

**HEADWORK SEALING:**

- Brass

**MATERIALS:**

- CERAMIC temperature working range: 0-70°

**CONNECTION TAP:**


- Flexibles F3/8"

**STANDARD AND NORMS:**

- EN13792:2000
- DIN12898
- DIN12918
- ISO228/1

**OPTIONS:**

- With chemical nickel plating 17 microns inside



**WEIGHT:**

- 2156 gr.

**WARRANTY 2 YEARS**



**Nuova FAR  
Rubinetterie**

Disegnato da

Davide Riva

Approvato da

Michele Albertone

Verificato da

Michele Albertone

Data di emissione: 16/02/2012

Peso:

Scala: 1:5

Il presente disegno a termine di legge non può essere riprodotto o reso pubblico senza la nostra autorizzazione scritta.

**DISEGNATO CON CAD NON MODIFICABILE A MANO**